


Data sheet Otocompact Professional



Technical details

NAME	VALUE
Power supply tension	220/240 V(110/115 VAC)
Frequency	50/60 Hz
Absorbed power (MAX maximum configuration)	1200 VA
Weight	180 Kg
Maximum number of instruments	4
Insufflator air capacity	26 l/m (variation up to 250 Kpa)
Washing stringe capacity	350/400 cc/min (2 Bar)
Nominal suction (oil motor)	-86 KPa (-650 mm Hg,50L/min)
Power supply fuses	F1, F2 = 6,3A T
Fuse external Plug	F3 =4A T

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it