


Data sheet OtoP/V Professional



Technical details

NAME	VALUE
Maximum seat height	750 mm
Base lenght	740 mm
Chair weight	80 Kg
Fuses protection	1,6 A
Absorbed power	90 V.A
Lifting Capacity	180 Kg
Backrest push	800 N
Lifting column push	2000 N
Low tension column	24 Vdc
Power supply tension	230 V-50/60 Hz

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it