

Data sheet Inventis mod. Flute: Inventis middle ear analyzer mod. Flute

Technical details - Sistemi di misura / Measurement systems

NAME	VALUE
Man pump contr.	available in Plus model
Range press	da -600 a +400 daPa
Range compl	0.2 a 8.0 ml
Toni sonda	226 Hz a 85 dB SPL - 1000 Hz a 75 dB SPL - (only mod. HF)

Technical details - Timpanometria / Tympanometry

NAME	VALUE
Class	Class IIa (MDD 93/42)
Standards	Tympan.: EN 60645-5 /ANSI 5339, Type 2 - Safety: EN 60601-1 - EMC: EN 60601-1-2
Weight	2.0 kg
Dimension	32 x 32 x 9 cm
Software compatibility	Inventis Daisy o Noah (with dedicated moduls)
PC interface	USB (without driver)
Printer (optional)	Integrated thermal
Resolution	4,3", 480 x 272
Screen	Color graphical
Memory	1 tympan. + 8 reflex + 4 decay record. per ear
Test ETF	Intact eardrum (Plus e HF) - Perforated eardrum (solo HF)
Automatic tests	Quick A e B (Seq. di temp. e rifl)
Stim. duration	da 0.5 a 2.5 sec selecting, Reflex decay rif. 10 o 20 sec
Max intensity	110 dB HL (ipsi); 120 dB HL (contra)
Stimuli	250 (C), 500, 1k, 2k, 4k Hz - Solo Plus: 3k, 6k, 8k, (C) Hz, LPN, HPN, BBN, no stim.
Available test	Autothreshold, Fixed intensity, Growing intensity-Only in Plus model: Manual test, Reflex decay Soglia autom.
Stimulation	Ipsi e controlaterale

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it

Data sheet Inventis mod. Flute: Inventis middle ear analyzer mod. Flute

Press change	50, 100, 200, 300 daPa/sec, AUTO
Type of test	Automatic and manual tympanometry

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
Capitale sociale: 50.000,00 Rea: BO-515493
Partiva Iva e CF: 03390661209
Pec: medicaresolutions@legalmail.it