


## Data sheet OtoP/V Elegance



### Technical details

| NAME                 | VALUE          |
|----------------------|----------------|
| Power supply tension | 230 V-50/60 Hz |
| Low tension column   | 24 Vdc         |
| Lifting column push  | 2000 N         |
| Backrest push        | 800 N          |
| Lifting Capacity     | 145 Kg         |
| Absorbed power       | 90 V.A.        |
| Fuses protection     | 1,6 A          |
| Chair weight         | 80 Kg          |
| Base leght           | 740 mm         |
| Maximum seat height  | 750 mm         |

### MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy  
 Capitale sociale: 50.000,00 Rea: BO-515493  
 Partiva Iva e CF: 03390661209  
 Pec: medicaresolutions@legalmail.it