


Data sheet OtoP/i Elegance



Technical details

NAME	VALUE
Power supply tension	NC - hydraulic piston
Hydraulic pump pus	3000 N
Backrest push	800 N
Lifting Capacity	150 Kg.
Chair weight	80 Kg
Base lenght	740 mm
Maximum legrest height	680 mm

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it