


Data sheet Otocompact Steel



Technical details

NAME	VALUE
Power supply tension	230 V +/- 10% optional 110V +/-10%
Frequency	50/60 Hz
Absorbed power (MAX maximum configuration)	1200 VA
Weight	120 Kg
Maximum number of instruments	4
Insufflator air capacity	26 l/m (variance up to 250 Kpa)
Washing stringe capacity	350/400 cc/min (2,2 Bar)
Nominal suction (oil motor)	-86 KPa (-650 mmHg) 50 L/min with infinite variation
Power supply fuses	F1, F2 = 6,3 AT
External plug fuse	F3 = 4AT

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it