


Data sheet Otoplus



Technical details

NAME	VALUE
Power supply tension	220/240 V (110/115 VAC)
Frequency	50/60 Hz
Absorbed power (MAX configuration)	1200 VA
Weight	100 Kg
Maximum number of instruments nr.	4
Insufflator air capacity	26 l/m (variance up to 250 Kpa)
Washing stringe capacity	350 cc/min
Nominal suction (oil motor)	-82 Kpa (-615mm Hg, 34lt/min)
Power supply fuses	F1, F2 = 6,3A T
Fuse external Plug	F3 =4A T

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it