


Data sheet Steel TS-X



Technical details

NAME	VALUE
Feeding tension	115/230 V
Frequency	50/60 Hz
Power absorbed (according to configuration)	650 VA
Weight	100 Kg
Maximum number of instruments	4
Maximum air pressure	3.5 Bar
Air capacity	3,5 l/m
Maximum water pressure	2.2 Bar
Feeding Fuse	F1,F2 = 6,3A T
24 Vdc fuses	F3=5 AT

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it