

Data sheet Sterix 18

Technical details

NAME	VALUE
Chamber capacity	18 lt
Chamber size	Ø 260 x D 365 mm
External Dimensions	L 495 x H 545 X D 540 mm
Min Top Depth	50 cm
Supply Voltage	230V 50/60Hz
Max Power Consumption	2.800 W
Next Weight	59
Heating System	Adaptive Integrated Heat + double steam generator
Vacuum pump	Tween Head 44lt/m > 0,97 bar
Reservoirs	2 inspectionable 5 lt tanks
Water quality control	Yes
Drain front/rear	Yes
WI-FI Module	Yes
IOT (Internet of things)	Yes
Printer connection	Yes
NFC	Yes
Night cycle	Yes
Warranty	2 years or 1500 cycles - 20 years for the chamber

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it