

Data sheet Stroboscopic light sources

Technical details - Basic version

NAME	VALUE
Light source	with LED on handpiece
Color temp	6,500°K
LED intensity	800 lumens
LED lifetime	50,000 hours
Mode of use	continuous and stroboscopic light
F0 ex-traction	manual or automatic
F0 range	70 Hz - 1,000 Hz
Slow motion	3 speeds
Microphone	clip-on and contact microphones included
Footswitch	Single pedal footswitch included
Light cable conn	Storz compatible, other standards available on request
Outputs	audio and USB
Weight	5 Kg / 11 lbs
Dimensions	180x100x140 mm (WxHxD)
Standards	Safety: EN 60601-1EMC: EN 60601-1-2

Technical details - Plus version

NAME	VALUE
Standards	Safety: EN 60601-1EMC: EN 60601-1-2
Dimensions	300x110x290 mm (WxHxD)
Weight	5 Kg / 11 lbs
Light cable conn	Storz compatible, other standards available on request
Footswitch	2 pedals footswitch included
Microphone	clip-on and contact microphones included
Slow motion	3 speeds
F0 range	70 Hz - 1,100 Hz

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
 Capitale sociale: 50.000,00 Rea: BO-515493
 Partiva Iva e CF: 03390661209
 Pec: medicaresolutions@legalmail.it

Data sheet Stroboscopic light sources

F0 ex-traction	manual or automatic
Mode of use	continuous and stroboscopic light
LED lifetime	20,000 hours
LED intensity	2,000 lumens
Color temp	6,800°K
Light source	white LED, embedded in the device housing

MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy
Capitale sociale: 50.000,00 Rea: BO-515493
Partiva Iva e CF: 03390661209
Pec: medicaresolutions@legalmail.it