


## Data sheet Suprema U-500



### Technical details

| NAME                                         | VALUE                |
|----------------------------------------------|----------------------|
| Power supply tension                         | 100-240 VAC 50/60 Hz |
| Fuses                                        | F1=F2=1,6 AT         |
| Max absorbed power                           | 200VA                |
| Class according to IEC 60601-1               | IB                   |
| Sectionable power feeding cable              | L=3mt                |
| Lifting motor                                | 8000N                |
| Backrest motor                               | 6000N                |
| Trendelenburg motor                          | 6000N                |
| Lifting capacity                             | 200 Kg               |
| Chair weight                                 | 180 Kg               |
| Lifting excursion, Trendelenburg inclination | 350 mm, 20° +/- 2°   |

### MEDI-CARE SOLUTIONS S.R.L

Via Pietro Nenni, 3 - 40026 Imola (BO) - Italy  
 Capitale sociale: 50.000,00 Rea: BO-515493  
 Partiva Iva e CF: 03390661209  
 Pec: medicaresolutions@legalmail.it